


Sizing Guide: Oxford® Universal Loop Sling MultiFit

	Height (cm)	Girth (cm)	Max SWL	Size	Colour Code
Paediatric (XS - S)	75 - 95	45 - 55	227 kg / 35st	1 (XS)	
	95 - 110	50 - 60	227 kg / 35st	2 (XS)	
	110 - 130	55 - 70	227 kg / 35st	3 (S)	
	130 - 150	65 - 80	227 kg / 35st	4 (S)	
Adult Male (M - XL)	135 - 160	70 - 85	227 kg / 35st	5 (M)	
	145 - 170	75 - 90	227 kg / 35st	6 (M)	
	155 - 180	85 - 100	227 kg / 35st	7 (M)	
	160 - 190	95 - 115	227 kg / 35st	8 (M)	
	160 - 190	105 - 130	227 kg / 35st	9 (L)	
	160 - 190	125 - 150	227 kg / 35st	10 (L)	
	160 - 190	145 - 175	284 kg / 44st	11 (XL)	
	160 - 190	170 - 200	284 kg / 44st	12 (XL)	
Adult Female (M - XL)	130 - 155	85 - 100	227 kg / 35st	5 (M)	
	135 - 160	95 - 110	227 kg / 35st	6 (M)	
	140 - 165	105 - 125	227 kg / 35st	7 (M)	
	140 - 165	120 - 140	227 kg / 35st	8 (M)	
	140 - 165	135 - 155	227 kg / 35st	9 (L)	
	140 - 165	150 - 175	227 kg / 35st	10 (L)	
	140 - 165	170 - 200	284 kg / 44st	11 (XL)	



Measuring the Client

You will need to take two measurements to find the most suitable sling size for your client; height and girth. See guidance below whilst referring to the diagram and table above for assistance.

1. Height

Take your client's full height as shown. It may prove more accurate to measure height from a recumbent (lying) position.

2. Girth

For female clients, measure around the hips. For male clients and children, measure around the chest just under the arms.

Universal MultiFit

This Universal MultiFit incorporates over-sized leg supports to specifically cater for clients with larger than average hips/thighs. Should you find a client's girth size is too large to match the height in the table above, try these for comparison.